

UNITED STATES TRANSURANIUM AND URANIUM REGISTRIES  
ANALYTICAL PROCEDURE MANUAL

**USTUR 200: ANION EXCHANGE ISOLATION OF PLUTONIUM FROM PREPARED TISSUE SOLUTIONS**

<b>Purpose</b>	Anion exchange for isotopes of plutonium	<b>Method Number</b>	USTUR 200
<b>Original Date</b>	10/10/95	<b>Author</b>	Radiochemistry Staff
<b>Revision Number</b>	3	<b>Approved By</b>	James T. Elliston
<b>Revision Date</b>	3/1/00	<b>Approval Date</b>	1/31/01

**SAFETY NOTE:** Before beginning this procedure, read all of the Material Safety Data Sheets for the chemicals listed in Section 5 of this procedure.

**1. Principle of Method**

- 1.1. Radionuclides of plutonium are isolated from ashed and acid-dissolved tissue samples (USTUR 100).
- 1.2. An aliquot of the tissue sample solution is selected for analysis on the basis of wet sample mass, ash weight, and estimated sample alpha activity.  
  

**NOTE:** A concentration of 25 mg ash per mL of 9 M HCl-Cl<sub>2</sub> should not be exceeded. Ash weights for bones are recorded, while ash weights of soft tissues are assumed to be 1% of wet weights.
- 1.3. Plutonium-242 tracer is added to the aliquot. Plutonium is isolated by anion exchange and is electrodeposited onto stainless-steel disks.
- 1.4. Concentrations of <sup>238</sup>Pu and <sup>239</sup>Pu are determined by alpha spectrometry.
- 1.5. Chemical losses are identified and corrected on the basis of the recovery of the added <sup>242</sup>Pu tracer.

**NOTE:** If only americium is being analyzed, use procedure USTUR 210.

**2. Minimum Detectable Activity (MDA)**

- 2.1. MDA is limited by the counter background.
- 2.2. For routine measurements with a 75,000-s sample counting period; 600,000-s background count (20 counts), 60% tracer recovery, and a counting efficiency of 25%, the L<sub>D</sub> is 0.001 Bq per sample (0.05 dis/min) of <sup>238,239</sup>Pu. See USTUR 600 for alpha spectrometry procedures.

### 3. Accuracy and Precision

- 3.1. Average blank recovery of  $^{242}\text{Pu}$  tracer is equal to  $87\% \pm 8\%$  at the 0.15-Bq (9 dpm) concentration level as of 6/10/96.

### 4. Apparatus

- 4.1. Hot plates.
- 4.2. Magnetic stirring hot plates.
- 4.3. Stirring bars: Teflon-coated.
- 4.4. Hot plate thermometer with a range to  $200^{\circ}$ .
- 4.5. Watch glasses: assorted sizes to fit beakers used.
- 4.6. Bio-Rad Ion-exchange columns (Fig. 1): borosilicate glass barrel with polypropylene reservoir, column tip, and bed support, 20 cm long by 1.0 cm i.d.
- 4.7. Rack: to support ion-exchange columns.
- 4.8. Glass beads: 2 or 3 mm diameter.
- 4.9. Graduated cylinders: various sizes.
- 4.10. Beakers: various sizes.
- 4.11. Wash bottle, 500 mL.

### 5. Reagents

- 5.1. 18 M  $\Omega$  (nanopure) deionized water
- 5.2. Hydrochloric acid (concentrated, 36.5 – 38%, reagent grade).
- 5.3. Hydrogen peroxide (30%, reagent grade).
- 5.4. Chlorinated hydrochloric acid solution (9 M HCl-Cl<sub>2</sub>). Add 5 mL H<sub>2</sub>O<sub>2</sub> (30%) and 750 mL concentrated HCl to 250 mL nanopure H<sub>2</sub>O. Heat for 1 hour at  $120^{\circ}\text{C}$ , noting formation of yellow color. Make fresh just before use.
- 5.5. Hydrochloric acid (0.6 M). Add 50 mL of concentrated HCl to 950 mL of nanopure H<sub>2</sub>O.

UNITED STATES TRANSURANIUM AND URANIUM REGISTRIES  
ANALYTICAL PROCEDURE MANUAL

- 5.6. Bio-Rad anion-exchange resin (AG 1-X4, 100-200 mesh) chloride form. Make a slurry of half resin, half nanopure water in a wash bottle.
- 5.7. Nitric acid (concentrated, 69 – 71%, reagent grade).
- 5.8. Nitric acid (8M). Add concentrated HNO<sub>3</sub> to an equal volume of nanopure H<sub>2</sub>O.
- 5.9. Hydrochloric acid (1.2 M). Add 100 mL of concentrated HCl to 900 mL nanopure H<sub>2</sub>O.

## 6. Procedure

- 6.1. Sample preparation.
  - 6.1.1. Refer to procedure USTUR 100, Tissue Ashing, Sample Dissolution, Sample Aliquot Selection, and Tracer Addition.
  - 6.1.2. Do not add U-232 tracer to an aliquot to be analyzed for Pu. The order of stripping will allow Th-228 (U-232 daughter) to interfere with Pu-238. An alternate elution step must be used!
- 6.2. Sample dissolution.
  - 6.2.1. To the previously prepared sample, add enough 9 M HCl-Cl<sub>2</sub> (a minimum of 100 mL is needed, keeping a concentration of 25 mg ash per mL or less) to completely dissolve the sample when heated to 120°C on a hot plate. Use stirring hot plates with magnetic stir bars to hasten the dissolution of large samples. Cover the beaker with a plain watch glass while heating.
  - 6.2.2. When the sample is dissolved, heat it for an additional 5-10 min on a hot plate with a surface temperature of 120°C. Remove the sample from the hot plate and allow it to cool to room temperature.
- 6.3. Anion exchange separation.
  - 6.3.1. Add 12-mL of nanopure water to the column and mark the column at the water level. Drain off the H<sub>2</sub>O.
  - 6.3.2. Fill the column with a slurry of AG 1-X4 resin to a settled volume of slightly more than 12 mL (approx. 0.5 cm above the mark).
  - 6.3.3. Rinse all of the resin down from the sides of the column using nanopure water, then 0.6 M HCl. Once all of the acid has drained, add approximately 1 cm depth of glass beads on top of the resin to prevent

UNITED STATES TRANSURANIUM AND URANIUM REGISTRIES  
ANALYTICAL PROCEDURE MANUAL

disturbance as reagents are added. Place a reservoir on the top of the column.

- 6.3.4. Wash the column with 60 mL (5 column volumes) of 9 M HCl-Cl<sub>2</sub>, allowing the solution to drip into a waste beaker. Discard this solution into a hazardous waste container (Non-rad).
- 6.3.5. If americium is also to be determined in the sample, place a clean labelled collection beaker under the column. This must be large enough to contain the sample volume plus the 85 mL (7 column volumes) of 9 M HCl-Cl<sub>2</sub> wash solution that will follow.
- 6.3.6. Transfer the sample to the ion exchange column. Rinse the beaker twice with 9 M HCl-Cl<sub>2</sub>, adding the washings to the column.
- 6.3.7. Adjust the flow rate to between 1 and 3 mL/min.
- 6.3.8. When the solution has drained to the top of the resin bed, add 85 mL (7 column volumes) of 9 M HCl-Cl<sub>2</sub> to the column.
- 6.3.9. When the solution has reached the top of the resin bed, remove the americium collection beaker, if used, and replace it with a waste beaker for the next wash. Evaporate the americium fraction to dryness on a hot plate with a surface temperature not greater than 120°C.

**NOTE:** Proceed to procedure USTUR 300, "Anion Exchange Isolation of Americium from Prepared Tissue Solutions," for the determination of americium.

- 6.3.10. Wash the column with 120 mL (10 column volumes) of 8 M HNO<sub>3</sub>. If the sample is also to be analyzed for uranium, discard the first 30 mL (2.5 column volumes) of this wash into a radioactive waste container. Replace the waste beaker with a clean labeled 150-mL beaker and collect the remaining 90 mL (7.5 column volumes) for the uranium procedure, USTUR 400. Evaporate this uranium fraction to dryness on a hot plate set to 150°C, but do not bake.

**NOTE:** Washing with 8 M HNO<sub>3</sub> removes the iron from the resin bed. At this point the deep yellow color is removed from the column.

- 6.3.11. When the wash solution has reached the top of the resin bed, replace the waste beaker with a labeled 150-mL beaker for the collection of the plutonium labeled fraction.
- 6.3.12. Elute the plutonium with 85 mL (7 column volumes) of 1.2 M HCl containing 2 mL 30% H<sub>2</sub>O<sub>2</sub> per 100 mL 1.2 M HCl.

- 6.3.13. Evaporate the eluate to dryness on a hot plate set to 120°C.
- 6.3.14. Add 5-10 mL of concentrated HNO<sub>3</sub> to the beaker and evaporate to dryness at 150°C. Repeat this step as necessary until no (or little) residue is present.
- 6.4. Proceed to procedure USTUR 500, "Electrodeposition of Americium, Plutonium, and Uranium" or USTUR 510.

## 7. Source Materials

- 7.1. I.K. Kressin, *Anal. Chem.* **49**, 842 (1977).
- 7.2. I.K. Kressin, W.D. Moss, E.E. Campbell, and H.F. Schulte, *Health Phys.* **28**, 41 (1975).
- 7.3. LANL method RT 400.

UNITED STATES TRANSURANIUM AND URANIUM REGISTRIES  
ANALYTICAL PROCEDURE MANUAL




Figure 1. Ion Exchange Column