

UNITED STATES TRANSURANIUM AND URANIUM REGISTRIES
ANALYTICAL PROCEDURE MANUAL

**USTUR 220: SEPARATION AND PURIFICATION OF PLUTONIUM FROM
AMERICIUM IN PRE-CONCENTRATED SAMPLES**

Purpose	Separation and purification of plutonium from Americium in pre-concentrated samples	Method Number	USTUR 220
Original Date	11/15/99	Author	Dorothy Stuit
Revision Number	2	Approved By	Jim Elliston
Revision Date	05/12/03	Approval Date	05/30/03

SAFETY NOTE: Before beginning this procedure, read all of the Material Safety Data Sheets for the chemicals listed in Section 3 of this procedure.

1. Principle of Method

- 1.1. Plutonium and americium (pre-concentrated from tissue samples in USTUR 150) are quantitatively separated by anion exchange.
- 1.2. The plutonium portion is prepared for electrodeposition (USTUR 510).
- 1.3. The americium portion is purified through further separation in procedure USTUR 310.

2. Minimum Detectable Activity (MDA)

Not applicable.

3. Accuracy and Precision

Not applicable.

4. Apparatus

- 4.1. Hot plates.
- 4.2. Watch glasses: assorted sizes to fit beakers used.
- 4.3. Hot plate thermometer with range to 200°C.
- 4.4. Bio-Rad Ion-exchange columns (Fig. 1): borosilicate glass barrel with polypropylene reservoir, column tip, stop-cock and bed support; 20 cm long by 1.0 cm i.d. or equivalent.

UNITED STATES TRANSURANIUM AND URANIUM REGISTRIES
ANALYTICAL PROCEDURE MANUAL

- 4.5. Rack: to support ion-exchange columns.
- 4.6. Glass beads: 2 or 3 mm diameter.
- 4.7. Graduated cylinders.
- 4.8. Beakers: various sizes.
- 4.9. Metric ruler.
- 4.10. Wash bottles: 500 mL.

5. Reagents

- 5.1. 18 M Ω deionized water (D.I. water).
- 5.2. Nitric acid (concentrated 69-71%, reagent-grade).
- 5.3. Nitric acid (8M). Dilute concentrated HNO_3 with an equal volume of D.I. water by adding the acid to the water.
- 5.4. Sodium nitrite (reagent-grade).
- 5.5. Bio-Rad anion exchange resin (AG1-X4, 100-200 mesh) chloride form. Make a slurry of half resin, half D.I. water in a wash bottle.
- 5.6. Hydrochloric acid (concentrated 36.5-38%, reagent-grade).
- 5.7. Hydrochloric acid (9M). Add 750 mL of concentrated HCl to 250 mL of D.I. water.
- 5.8. Hydrochloric acid (0.6 M). Add 50 mL of concentrated hydrochloric acid to 950 mL of D.I. water.
- 5.9. Hydroxylamine hydrochloride (reagent grade crystals).
- 5.10. Ammonium iodide (reagent grade crystals).
- 5.11. 0.6 M hydrochloric acid/0.05 M ammonium iodide. Dissolve 0.72 g of NH_4I in 0.6 M HCl. Dilute to 100 mL with 0.6 M HCl. Make fresh just prior to use.

6. Procedure

- 6.1. Sample preparation.

UNITED STATES TRANSURANIUM AND URANIUM REGISTRIES
ANALYTICAL PROCEDURE MANUAL

6.1.1 Sample aliquot is obtained from USTUR 150, step 6.4.15.

6.2 Sample dissolution.

6.2.1. Add 50 mL of 8 M HNO₃ to sample aliquot and cover with watchglass. Warm on 150°C hotplate, if necessary, to dissolve sample. Remove samples from the hotplates.

6.2.2. Turn hotplate up to 300-350°C (with samples off hotplate).

6.2.3. Add approximately 0.1 g sodium nitrite (tip of large spatula full) to each sample. Swirl sample, replace watchglass, and place on hotplate. Heat until just starting to boil.

NOTE: Solution should be maintained at a slight yellow color. Do not heat until clear! Once sodium nitrite has been added to the samples, the separation must be completed in a timely manner. Columns should not be left dry for longer than 20-30 minutes.

6.2.4. Remove sample from hotplate and bring to room temperature before proceeding with anion exchange.

6.3. Anion exchange separation.

6.3.1. Using a ruler, mark the ion exchange column 8 cm above the frit.

6.3.2. Fill the column with the slurry of AG1-X4 resin to a settled depth at the 8g cm mark.

6.3.3. Rinse all the resin down from the sides of the column using deionized water, then 8 M HNO₃. Once all of the acid has drained, add approximately 1 cm depth of glass beads on top of the resin.

6.3.4. Wash the column with 50 mL of 8 M HNO₃. Discard the wash in a hazardous waste acid container.

6.3.5. Place a clean 250 mL beaker labeled with sample number and "Am" to catch the americium portion of the sample.

6.3.6. Add cooled sample to column reservoir. Rinse beaker twice with approximately 5 mL 8 M HNO₃, adding the rinse to the column. Allow column to drain completely.

6.3.7. Rinse the column with 60 mL of 8 M HNO₃, collecting the column effluent in the beaker for americium analysis. Allow column to drain completely.

- 6.3.8. Remove americium beaker and bring to dryness on a hotplate set at 120°C. Cover with a watch glass or parafilm and store for purification by USTUR 310.
- 6.3.9. With a waste beaker under the column, rinse column with 30 mL of 9 M HCl. Allow column to drain completely. Discard effluent in a Radioactive waste container.
- 6.3.10. Place a clean 150 mL beaker labeled with sample number and “Pu” under column to catch plutonium portion of sample.
- 6.3.11. Cover tops of glass beads (in column) with hydroxylamine hydrochloride crystals.
- 6.3.12. Add 30 mL of 0.6 M HCl to the column and allow it to drain completely.
- 6.3.13. Add 20 mL of 0.6 M HCl/0.05 M NH₄I solution to the column and allow it to drain completely.
- 6.3.14. Add 20 mL of concentrated HNO₃ to plutonium beaker and bring solution to dryness on a hotplate set at 120°C.
- 6.3.15. Continue with step 5.1 in USTUR 510 to prepare the plutonium portion for electrodeposition.
- 6.3.16. Rinse the column with 50 mL of D.I. water. Discard the rinse into a radioactive waste container. Place the used resin into a resin collection beaker.

7. Source Materials

- 7.1. “Isolation of Plutonium, Americium, and Uranium from Urine, by Precipitation, Ion Exchange, and Extraction,” 4-16200-RHL-0037, Rocky Flats Plant.

UNITED STATES TRANSURANIUM AND URANIUM REGISTRIES
ANALYTICAL PROCEDURE MANUAL




Figure 1. Ion Exchange Column