

EdU detection in Formalin-fixed paraffin sections with ClickiT AlexaFluor 488 kit

1. Deparaffinize:

Xylene 2 X 5 min
100% EtOH 2 X 3 min
95% EtOH 2 X 3 min
70% EtOH 3 min

Prepare

PBS+3%BSA

2. Rinse with water to remove EtOH.
3. Draw circle with pap pen
4. Block for 10 minutes with PBS+3% BSA
5. Prepare **Click-iT reaction cocktail** (need approx 100 µl/section);
use within 15 minutes of preparing it.
6. Per **1 ml**:
 - (1) **860 µl 1x Click iT reaction buffer (Component D)**
(diluted fresh with dH2O from 10x reaction buffer concentrate "Component D")
 - (2) **40 µl component E** (CuSO₄)
 - (3) **5 µl AlexaFluor 488 azide** (one that comes with kit; made by mixing 70 µl
component C (DMSO) with component B until dissolved)
(A10266; purchased separately MW = 861.04;
0.5 mg dissolved in 100 µl DMSO=5.81 mM stock; , store at -20 for up to 1 yr)
2 µl per 1 ml cocktail=11.6 µM final)
 - (4) **100 µl Reaction buffer additive (Component F)** (diluted fresh from
10x stock with dH2O;
10x stock is made by adding 2 ml dH2O to Component F vial,
mixing until dissolved; store at -20 for up to 1 yr.
Once 10x stock is diluted, reaction buffer additive must be used that day)
7. Add **Click-iT reaction cocktail** to sections; incubate for 30 minutes at room temp,
protected from light
8. Wash 1x for 5 min with PBS+3%BSA
9. Wash 1x for 5 min with PBS
10. Dilute Hoescht 33342 (**Component G**) 1:2000 in PBS, apply to sections for 30 min,
protected from light. Remove Hoescht
11. Thaw ProLong Gold (No DAPI)
12. Wash 3x for 5 min with PBS
13. Mount with ProLong Gold (no dapi), cure in dark for 24h
Use this instead of permount