

“Low Level Sulfur Detection”

Quantitation of Low-Level Sulfur in Middle Distillate Fuels by ICP-MS Project 105

Lead investigator: Joe Gord, University of Dayton Research Institute
Project manager: Ana Gabrielian, FAA

October 13, 2025
Alexandria, VA

This research was funded by the U.S. Federal Aviation Administration Office of Environment and Energy through ASCENT, the FAA Center of Excellence for Alternative Jet Fuels and the Environment, Project 105 through FAA Award Number 13-C-AJFE-UD under the supervision of Ana Gabrielian. Any opinions, findings, conclusions or recommendations expressed in this material are those of the authors and do not necessarily reflect the views of the FAA.



Project 105

Low-Level Sulfur Detection

University of Dayton Research Institute

PI: Zach West, Joe Gord

PM: Ana Gabrielian

Cost Share Partner(s): Airlines, OEMs, Fuel Producers



Objective:

Development of an ICP-MS method capable of quantitating low levels of sulfur in middle distillate fuels—Target concentration range from 10s of parts-per-billion (ppb) to 10s of parts-per-million (ppm)

Project Benefits:

Improved sulfur detection limits over currently available methods supporting SATF candidate fuel testing and approval, advancement of low-sulfur fuel standards, and generation of data supporting research efforts such as CLEEN with partners like NASA and Boeing.

Durability, Performance, Health, and Stealth

Research Approach:

- Measuring sulfur in fuel samples by ICP-MS is especially difficult for several reasons including:
 - Low ionization efficiency for sulfur in an argon plasma
 - High coking potential for fuel samples
 - Significant isobaric interference
- Collision Cell Technology (CCT) is implemented to react sulfur ions with oxygen gas to form sulfur oxide ions which can be detected in a less-crowded mass channel ($m/z = 48$)
- Using ASTM D8110-25 as a guide for developing an ASTM-Approvable ICP-MS method.

Major Accomplishments (to date):

Draft method established with good performance across a variety of fuels in the 100s of ppb to low-ppm range.

Future Work / Schedule:

- Refine method, finalize performance statistics, conduct multi-lab robustness study, apply method to broader range of SATF/SBC fuels
- Incorporate non-interfering internal standard
- Publish as official UDRI Method


Introduction – why are we doing this?

- Reduced sulfur content in jet fuel supports durability, performance, health, and stealth.
- Legacy Sulfur Test Methods and Specifications only reach a few ppm (mg/kg).
- Existing ICP methods don't measure sulfur.
- Think of this as a research tool.

Method Number	Lower Limit [S]	Measurement	Use Case (Sulfur)
D5453	1 mg/kg	UV Fluorescence	D1655, D7566
D2622	3 mg/kg	Wavelength Dispersive X-Ray Fluorescence (XRF)	D1655, D7566
D4294	16 mg/kg	Energy Dispersive XRF	D1655, D7566
D1266	100 (5*) mg/kg	Lamp	Withdrawn 2025
IP 336	300 mg/kg	Energy Dispersive XRF	D1655, D7566
D7111	N/A	ICP-Atomic Emission Spectroscopy	N/A – Metals
D8110	N/A	ICP-Mass Spec	N/A – Metals



Introduction: Technique Overview




- Sulfur has high ionization energy (10.4 eV)
 - Only ~10% ionized in “standard plasma”
- Sulfur is light and inefficiently transported by ion optics.
- High-Carbon Matrix (Jet Fuel) requires addition of $^{16}\text{O}_2$ ($m/z = 32$)
 - $^{16}\text{O}_2$ = Isobaric Interference with ^{32}S
- $^{32}\text{S}^{16}\text{O}^+$ ($m/z = 48$) easy to make
 - Not obscured by $^{16}\text{O}_2$
 - $^{36}\text{Ar} + ^{12}\text{C}$ ($m/z 48$) *may* interfere

4




Show me the Data!

- Developed ICP-MS method to detect S as SO
- 6 Cal Standards (100 – 3000 ppb)
- Good precision (repeatability)
- Accuracy & Robustness testing underway



Calibration standards are made by diluting a Specpure 5000 $\mu\text{g/g}$ Sulfur (Dibutyl Sulfide) standard in Kerosene.




*Calculated by instrument software package




Sulfur Calibration – Spike Test


Accuracy




Precision



Application Across a Range of Test Fuels



EcoDemonstrator Test Case



- HEFA-SPK and Jet-A data points are all the same fuel with aliquots taken from different sampling locations.
- Fuels with D5453 data presented as 1000 ppb [S] were reported as <1000 ppb, i.e., below the quant limit.
- ICP-MS appears to consistently over-report compared to D5453.
- ICP-MS appears to differentiate between points that would not be differentiable based on D5453 Reproducibility data.

Y values are $\pm 2\sigma$ | X values are $\pm R$ (Reproducibility), DF = 2



Summary: Method Highlights and Next Steps

- Calibrated Range = 100 – 3000 ppb
- Limit of Detection = 15 ppb
- Limit of Quantitation = 300 ppb



FAA CENTER OF EXCELLENCE FOR ALTERNATIVE JET FUELS & ENVIRONMENT



Summary: Method Highlights and Next Steps

- Calibrated Range = 100 – 3000 ppb
 - Limit of Detection = 15 ppb
 - Limit of Quantitation = 300 ppb
-
- How did we do?
 - Do you have any fuel we can test?
 - Would you like to try this at your place?



Acknowledgements

- FAA AEE for sponsorship & oversight—Ana Gabrielian & Anna Oldani

Participants

- UDRI Fuel Science Group Members:
Amanda Arts, Willie Steinecker, & Zach West (PI)

For more information please contact:

Joe Gord, Ph.D.

joe.gord@udri.udayton.edu

Research Chemist
Fuel Science Group
Fuels & Energetic Materials Division



FAA CENTER OF EXCELLENCE FOR ALTERNATIVE JET FUELS & ENVIRONMENT

