

Rapid IR Fuel Screening Project 25

Lead investigator: Ronald K. Hanson, Stanford University
Project manager: Bahman Habibzadeh, FAA

October 14-16, 2025
Alexandria, Virginia

This research was funded by the U.S. Federal Aviation Administration Office of Environment and Energy through ASCENT, the FAA Center of Excellence for Alternative Jet Fuels and the Environment, project 25 through FAA Award Number 13-C-AJFE-SU-027 under the supervision of Bahman Habibzadeh. Any opinions, findings, conclusions or recommendations expressed in this material are those of the authors and do not necessarily reflect the views of the FAA.



Project 25

Rapid IR Fuel Screening

Stanford University

PI: Ronald K. Hanson

PM: Bahman Habibzadeh

Cost Share Partner(s): Stanford University



Objective:

Develop a **compact, low-volume prescreening tool** for the prediction of physical and chemical properties of sustainable aviation fuels (SAFs) using Fourier-Transform Infrared (FTIR) spectrometry and advanced statistical analysis methods.

- Project Benefits:**
- The FTIR prescreening approach will **make the SAF design and approval process less costly and more efficient**
 - This low-volume (<1 mL) method yields insights that are **complementary to other prescreening approaches** (e.g., GCxGC).

- Research Approach:**
- Develop regularized linear and non-linear models that **correlate the physical, chemical, and combustion properties of a fuel** (e.g., boiling point, heat of combustion, flash point, etc.) **with its vapor-phase, 2-15 µm FTIR absorption spectrum.**
 - Apply these models **to predict physical and chemical properties of next-generation SAFs** and fuel components.

- Major Accomplishments (to date):**
- **Measured spectra of 38 World Jet Fuel Survey (WJFS) samples and predicted molecular weights**
 - Including WJFS samples in training set **improves model performance for similar fuels**
 - Implemented **facility improvements** and began preliminary **measurements of liquid-phase spectra**

- Future Work / Schedule:**
- Begin preliminary **measurements for combined liquid- and vapor-phase training set**
 - **Determine effect of phase** on the performance of property prediction models
 - Develop a **combined liquid-vapor model**


Introduction

Concept

- The vapor-phase IR absorption spectrum of a hydrocarbon fuel contains quantitative information about molecular structure and functional groups
- Statistical models can be used to infer the physical and chemical properties of fuels from this spectral information

Objective

Develop a compact, low-volume prescreening tool for the prediction of physical and chemical properties of sustainable aviation fuels (SAFs) using Fourier-Transform Infrared (FTIR) spectrometry and advanced statistical analysis methods




- Fuel structure is evident in the IR spectra; the location and magnitude of absorption features reflect the type and number of functional groups, and can be correlated with physical/chemical properties



Broad categories of spectral analysis strategies

- Spectrum-Property Correlation Models – Strategies 1 (linear) & 3 (non-linear)
- Spectrum-Composition Models – Strategy 2 (Functional groups)
Strategy 4 (Molecular species)




Past and Current Work

- ✓ Compared FTIR-based approaches with alternative methods in literature, revealing comparable performance
 - ✓ Successfully predicted the temperature-dependent thermodynamic properties for distillate fuels
 - ✓ Implemented facility improvements to allow for higher temperature and extended wavelength measurements
- *The following slides will focus on our recent progress analyzing model performance on 38 WJFS samples, training dataset populations, and preliminary comparison of vapor- and liquid-phase spectra*



FTIR Spectra of WJFS Samples




- 38 WJFS samples were measured in the vapor-phase
- Spectral variations prominent 3 and 7 μm regions
- Differences in spectra reveal functional group and compositional variation
 - Indicates opportunity for machine learning model to relate spectral variability to property variability

First application: Prediction of molecular weight...



Prediction of WJFS Molecular Weights




- **Initial training set** was developed for wide range of hydrocarbons and blends, but included only 8 conventional jet fuels
- Strategy 1 was trained on this limited set and applied to the 38 WJFS samples
- Overpredicts molecular weight with a Mean Absolute Error of 15.86 g/mol
- **Why:** Training set lacks sufficient number of samples with properties closely resembling conventional jet fuel!

Solution: Improve training set by incorporating more conventional fuels!



Analysis of Training Dataset Population




- WJFS samples show a relatively narrow distribution of molecular weight
- Initial training set spans a wider range and dips near the WJFS peak region
- Tailored set yields distribution that better represents conventional jet fuels

Emphasizes the need for a training set that reflects the target fuel types

Our broad dataset enables down-selection for tailored model development



Prediction of WJFS Molecular Weights



- 34 WJFS samples were added to initial set then the total set was **down-selected to molecular weight range of 135-185 g/mol**
- This **tailored training set** better aligns with conventional fuel range while preserving a robust sample size (total of 120 spectra)
- **4 WJFS samples kept out of training process** to simulate performance on future Jet-A samples
- Tailored model yields improved predictive performance with MAE of **0.9 g/mol** for **tailored training set** and **1.6 g/mol** for the **test set**

Size and distribution of the training set are critical for performance.

Performance improves with a larger, relevant training set.



Prediction of WJFS Molecular Weights

Conclusions

- Strategy 1 with limited initial training set overpredicts molecular weight of WJFS samples
- Adding WJFS samples and tailoring training set enhances model performance on conventional fuels
- Future work should expand and learn to tailor training sets for different fuel types

Next: Facility Improvements...




FAA CENTER OF EXCELLENCE FOR ALTERNATIVE JET FUELS & ENVIRONMENT



Facility Improvements



Heated Optical Cell



Nicolet iS50
FTIR Spectrometer




ATR Crystal

- Nicolet iS50 allows for expanded spectral range (1-200 μm), higher SNR, and better purge control and measurement of both vapor and liquid-phase
 - Heated optical cell allows for vapor-phase measurements at temperatures up to 300°C
 - ATR crystal cell used for measurement of liquid-phase samples



Comparison of Vapor and Liquid Spectra



Vapor-phase spectra

- Exhibit well-resolved fine features and follow straightforward mixing rules
- Useful for inferring composition

Liquid-phase spectra

- Exhibit broader, shifted bands
- Potentially useful for the inferring properties influenced by non-linear or intermolecular effects

Can we combine data for improved performance?



Comparison of Vapor and Liquid Spectra

Next Steps

- Compile a training set containing both vapor and liquid-phase spectra for each sample
- Determine effect of phase on the performance of property prediction models
- Develop a combined liquid-vapor model to achieve optimal performance for all properties



Acknowledgements

- Dr. Bahman Habibzadeh – Project management

Participants

- Prof. Ronald K. Hanson – Principal Investigator (rkhanson@stanford.edu)
- Dr. Jesse Streicher – Research Scientist (jessewst@stanford.edu)
- Dr. Miguel Figueroa Labastida – Postdoctoral Fellow (mfiglab@stanford.edu)
- Andrew Klingberg – Graduate Student (ark7@stanford.edu)



FAA CENTER OF EXCELLENCE FOR ALTERNATIVE JET FUELS & ENVIRONMENT

